

Picasso-Trio - 3-D Dental Imaging

Digital Panoramic & Cephalometric, CBCT X-ray Imaging System – 3 in 1 System

Model

Picasso-Trio - NP Panoramic only [CT upgradeable]

Picasso-Trio - NC Panoramic & Cephalometric [One sensor, CT upgradeable]
Picasso-Trio - SP8*5 Panoramic only with CT [Field of View (FOV) 8*5 cm]

Picasso-Trio - SC8*5 Panoramic & Cephalometric [One sensor] with CT [Field of View (FOV) 8*5 cm] - 3 in 1

Picasso-Trio - SP12*7 Panoramic only with CT [Field of View (FOV) 12*7 cm]

Picasso-Trio - SC12*7 Panoramic & Cephalometric [One sensor] with CT [Field of View (FOV) 12*7 cm] - 3 in 1

General Features

Picasso-Trio is the diagnostic system which consists of Digital Panoramic & Cephalometric Dental X-ray System, and Computed Tomography X-ray System with Cone Beam Technology. Especially, the advanced digital imaging process allows considerably efficient diagnosis, all kind of management information, real-time sharing of image and its information on network. It's equipped with the-state-of-the-art CT sensor to capture 3-D X-ray Computed Tomogram Scanned Image.

Movement Technology Multi-motor with digital trajectory control / motorized carriage movement

Patient Positioning

Panoramic Triple [Midsagittal / Frankfurt / Canine] laser beam positioning system
CT Triple [Midsagittal / Vertical / Horizontal] laser beam positioning system

Patient Positioning Aid Electric-motive chin rest, bite, mirror

Head-rest

Edentulous bite positioner

TMJ pointers

Ear rods for CEPH mode

Nasion support with vertical millimeter scale for CEPH mode

Voice Instruction Typical – English / Available Local Language


Up/Down Movement Silky up/down with one stage speed by electric-motor


Computed Tomography Examination Programs

- 3-D Volume Rendering
- Axial View
- Coronal View
- Sagittal View
- Cross-Sectional View
- Panoramic View
- MPR Multi-Planar Reconstruction
- Region of Interest [ROI] Analysis Bone Density [Hounsfield Units]
- Capturing Mode

: Mandibular(Mn), Occlusion(Occl), Maxillary(Mx), TMJ Right(Rx), TMJ Left(Tx) for FOV12*7







[Sample Image from FOV 12*7: Capturing Mode – Mandibular]

Standard Examination Programs

In each program the compensation for the spinal column is obtained by means of exposure parameter modulation, optimized in accordance with the selected anatomic program.

- Standard Panoramic
- Hemi-Panoramic (Left and Right)
- Frontal Dentition
- TMJ Open/Close mouth: 4 views [Right Open-Right Close-Left Close-Left Open] are taken on the same image
- Maxillary Sinus



Cephalometric Examination Programs

- Latero-Lateral
- Anterior-Posterior
- Carpus
- ServoMento Vertex







Anatomic Programs

CT Mode

• Patient: 3 choices: adult, weak, child

• Type: 2 choices: A-mode(normal), B-mode(sturdy short neck)

Panoramic & Cephalometric Mode

Patient Size:
 4 choices: adult man, adult woman, weak, child

Patient Type: 3 choices: hard, normal, soft
Arch Shape: 3 choices: wide, normal, narrow

Image Processing

CT Mode: Metal Artifact [Normal, Metal]
 Standard Examination Programs: Metal Artifact [Normal, Metal]

Image Reconstruction Time & File Size for CT Mode

• FOV 8*5

Voxel Size		File Size (MB)			
	Normal (scan time: 15 sec)		High (scan time: 24 sec)		
	Normal	Metal	Normal	Metal	
0.2	N/A	N/A	240	360	85
0.3	N/A	N/A	TBD	TBD	TBD

FOV 12*7

Voxel Size	Reconstruction Time (sec)				File Size (MB)
	Normal (scan time: 15 sec)		High (scan time: 24 sec)		
	Normal	Metal	Normal	Metal	
0.2	180	360	300	720	222
0.3	75	165	110	250	70

^{*} Image reconstruction time can be changed by computer specification and/or its working condition.

Image Magnification

• Computed Tomography Examination Programs 1.60 constant (virtual 1.00 constant)

Standard Examination Programs
 Cephalometric Examination Programs
 1.30 constant
 1.10 constant

FOD, ODD, FDD (mm)

Mode	FOD	ODD	FDD	
	(Focal Spot to Object Distance)	(Object to Detector Distance)	(Focal Spot to Detector Distance)	
CT	424.0	254.0	678.0	
Pano	422.7	173.0	595.7	
Ceph	1,580.7	221.0	1,801.7	


^{*} Metal Artifact: Please select "Metal" for patient with 2 ~ 3 metals, i.e. implants and/or crowns.



X-ray Generator

X-ray Beam Formation

Cone Beam



• High frequency generator, constant potential, micro processor controlled

Ripple < 4%

Inverter frequency
 36 kHz push-pull

• Tube type D-051, stationary anode type

Nominal power
 Below than 1.3 KW

Tube voltage

CT Mode 50 - 90 kV (adjustable by 1kV)Panorama Mode 40 - 90 kV (adjustable by 1kV)Cephalo Mode 40 - 90 kV (adjustable by 1kV)

Tube current

CT Mode 2.0 – 10.0 mA (adjustable by 0.1 mA)

Panorama Mode 2 – 10 mA (adjustable by 1 mA)
Cephalo Mode 2 – 10 mA (adjustable by 1 mA)

High voltage

DC

Exposure time

[Computed Tomography Examination Programs]

Scan Time

Normal Mode 15 sec (Available FOV 12*7 only)

High Mode 24 sec

[Panorama Examination Programs]

Standard Panoramic (Adult, Weak/Child) 13.2 sec/ 11.3 sec

Hemi-Panoramic (Right and Left) 6.6 sec Frontal Dentition 6.6 sec

TMJ Open/Close 14 sec (4 * 3.5 sec)

Maxillary Sinus 13.2 sec

[Cephalometric Examination Programs]

Latero-Lateral15 secAnterior-Posterior13.5 secCarpus13.5 secSMV13.5 sec

• Cooling by force, one minute for cooling / Protect > 50° C



X-ray Tube

Focal Spot

Heat storage capacity

Total filtration

0.5 mm

30 kJ (40 Khu)

2.8mm Al eq.

Collimator

Primary collimator

Fixed in basic configuration

Motorized positioning for PANO, CEPH & CT configuration

Secondary collimator at CEPH

Motorized positioning for CEPH configuration

Standard Accessories

- Bite blocks with supporters
- Disposable bag (Small & Large)



User Interface

All operating functions are easily controlled on the PC (kVp, mA, Image Capture Mode etc.) Voice instruction

[Computed Tomography Examinations Programs – FOV12*7]



Capturing Mode: Mandibular(Mn) – Occlusion(Occl) – Maxillary(Mx) – TMJ Right(Rx) – TMJ Left(Lx)

* NOTE: Above capture screen is subjected to FOV12*7.

[Standard Examination Programs - SC]



[Cephalo Examination Programs - SC]





Digital Image Acquisition System

Computed Tomography Image Detector

Technology
 High resolution flat panel detector

● Pixel size 200 µ

Voxel solution
 0.2 mm x 0.2 mm x 0.2 mm - Typical

available 0.3 mm3

Slices for 3-D image reconstruction
 306 (FOV12*7/Voxel 0.2 mm²)

206 (FOV12*7/Voxel 0.3 mm³)

Resolution of slice
 606 x 606 (FOV12*7/Voxel 0.2 mm²)

408 x 408 (FOV12*7/Voxel 0.3 mm³)

Projections

Normal Approx. 450

High Approx. 720

Field of view (FOV) 8*5 cm & 12*7 cm

Frame rate
 30 fps

Gray scale
 12 bits (Adjusted 16 bits for CT Number)

Slice width (gap)
 Sliced image thickness
 0.1 ~ 10 mm (default - 2 mm)
 0.1 ~ 10 mm (default - 1 mm)



Technology
 Digital multi-linear CCD sensor, safe detachable sensor with hand grip

Pixel size 48 μm

Active area

Panoramic mode 6.144 * 147.456 mm Cephalometric mode 6.144 * 221.184 mm


Image acquisition area

Panoramic mode 140.000 * 272.640 mm Cephalometric mode 201.216 * 216.000 mm

Gray level
 12 bits

Resolution

Panoramic mode 2840 * 1536 pixels Cephalometric mode 2096 * 2304 pixels



Axial Slices



EzImplant Basic [3-D Image Viewer]

EzImplant Basic is three-dimensional dental image viewer for prompt and accurate diagnosis with many useful functions as "various MPR function", "two-dimensional analysis" and "three-dimensional animating work" by loading DICOM format CT image, and more.

- Easy conversions through various rendering method as VR (Volume Rendering)/MIP/minIP/X-ray
- More accurate 3-D image by MPR rotating, curve, 3-D zoom rendering mode
- Cross-sectional view function for fast analysis
- Convenient management system of objects, color-map, opacity graph, preset files and more.
- Dicom format of reconstructed image for easy file saving and compacted interface that makes possible to operate with mouse button only in one window.
- Main Tool

View

Pan: Move image in the pane

Rotate: 3D/2D rotation

Zoom: Zoom in/out the image

Windowing: Adjust image brightness and coloring

Invert: Invert image brightness and coloring

Text overlay: Show image information

VOI (View of Interest) overlay

Measure

Distance: Ruler for 2 point, and Tapeline for various points

Angle

Profile: HU (Hounsfield Unit)

Area: Measure the area by drawing ROI (Region of Interest)

ROI (Region of Interest)

Segmentation

Draw Mask

3D Picker/Mask Overlay

Output

Capture: Pane, Region, Full screen

Print

CINE player: Moving clip

Task

MPR

Rotating axis: Move, rotate, adjust thickness

Oblique slice

Curve

Cross-sectional
General curve
Select latest item

Fine tuning

Opacity adjusting function

3D Edition and Magnification

Slice function, 3D distributing functions



EzImplant Professional [3-D Image Viewer]

EzImplant Professional supports following value-added functionality with all functionality at EzImplant Basic.

- Implant Simulation
- Canal Draw
- Show Bone Density
- Value-Added Utilities such as Free Draw, Report etc.

EasyDent [Viewer for Panoramic & Cephalometric Image]

- One click operation
- User friendly graphic interface
- Various image format

ВМР

Various image process & accurate measure tool

Patient and Image Archive

Main patient data

Examination for each patient

Radiological parameters set for each examination

Drawing

Various view modes

Implant simulation

Measure: 2 point distance, continuous distance, angle, user calibration

Image Processing

Bright adjustment

Contrast adjustment

Gamma adjustment

Coloring

Film coloring

Invert

Sharpen/Sharpen more

Local Enhancement (256x256, 128x128, 64x64)

Image rotation

Magnification etc.

Print

The system allows the user to print images on the screen with the following information:

Patient Information

Exposure parameters

Examination date

Images can be printed any Windows compatible printer, although image quality may differ as per the printer model.



System Requirements for the Computer

A. Image Capturing

Operating System
 Microsoft XP Professional Service Pack 2
 CPU
 Intel Pentium 4 950 3.4GHz Dual Core

Cash Memory
 4MB (2MB/Core)

• Chipset Intel 955X

HDD 250 GB (prefer bigger capacity)

Main Memory
 2 GB DDR2

Video Memory
 256 MB 128bit – Recommend NVIDIA Geforce

Network
Serial
2 EA of 10/100/1000 Gigabit
2 EA of RS232 serial port

DVD+/-RW

PCI Slot Min. 3EA - Mandatory
 Monitor Min, resolution 1280*1024

B. Image Viewer

Operating System
 Microsoft XP Professional Service Pack 2
 CPU
 Intel Pentium 4 950 3.4GHz Dual Core

Cash Memory
 Chipset
 Main Memory
 4MB (2MB/Core)
 Intel 955X
 1 GB DDR2

Video Memory
 256 MB 128bit - Recommend NVIDIA Geforce

Network
 1 EA of 10/100/1000 Gigabit

DVD+/-RW

Monitor
 Min, resolution 1280*1024

NOTE: E-WOO tested/verified HP Workstation XW4400 with CT system, model-named "Picasso - Trio", and we will not guarantee system performance by low specification of computer.



Mechanical Characteristics

Source to Image distance [Focal spot to Sensor]

CT 678.0 mm 591.1 mm Panoramic Cephalometric 1801.67 mm Max. 700.0 mm

Vertical column movement

Weight

Version without cephalometric unit 270.00 kg Version with cephalometric unit 310.00 kg


Max. 2312.0 mm Total height

Length * Width * Height

Version without cephalometric unit 1002(W) * 1476(D) * 2312(H) mm 2000(W) * 1476(D) * 2312(H) mm Version with cephalometric unit

Type of installation Base Stand

Dimensions (indicated in mm)





Electrical Characteristics

Power supply voltage
 AC 110/230V ± 10%

Frequency 50/60 HzPower rating 1.5KVA

Environmental Characteristics

Operating temperature 10 - 40°C
 Operating relative humidity 30 - 75%
 Operating atmospheric pressure 700 - 1060 hPa
 Transport and storage temperature -20 - 70°C


Transport and storage relative humidity
 < 90% non-condensing

Transport and storage atmospheric pressure
 500 – 1060 hPa

Standards and Regulations

This product is designed and produced to meet the following standards:

EN 60601-1, EN 60601-1-3, EN 60601-2-7, EN 60601-2-28, EN 60601-2-32, EN 60601-1-2, EN 61000-3-2, EN 61000-3-3, EN 61223-3-4, EN 61223-3-5, EN ISO 9001, EN ISO 13485



CE symbol grants the product compliance to the European Directive for Medical Devices 93/42 as a class IIB device. Authorized by **Grand-Duche De Luxemburg.**



Product Photography



NOTE: The actual specification may therefore be subjected to improvement and/or modification without notice.

E-WOO Technology Company Limited

138-2, Hagal-dong, Giheung-gu, Yongin-si, Gyeonggi-do, 446-930 Republic of Korea

Tel. +82-(0)31-288-1631 Fax. +82-(0)31-286-3787

E-mail. ct_technical@e-wootech.com URL. http://www.e-wootech.com